

Shipping

- Only the 50ml tubes with oral fluid need to be submitted to Genics ISO17025 Certified Laboratories
 - Samples to be submitted on the same day should be chilled and submitted with ice packs
 - Samples may also be frozen (at -20°C) and sent at your convenience
 - Submission form and instructions can be found at www.porkmultipath.com
- Samples should be send to:

Genics Head Office
Suite 14, 220 Boundary Street,
Spring Hill, QLD 4000 AUSTRALIA



Make Contact

If you have any questions or comments please contact:

info@genics.com
www.genics.com

Genics Head Office
Suite 14, 220 Boundary Street,
Spring Hill, QLD 4000 AUSTRALIA

© 2026 Genics Pty Ltd (Genics)
To the extent permitted by law, all rights are reserved, and no part of this publication covered by copyright may be reproduced or copied in any form or by any means except with the written permission of Genics. To the extent permitted by law, no part of this publication is to be interpreted by artificial intelligence.

Genics advises that the information contained in this publication comprises general statements based on scientific research. The reader is advised and needs to be aware that such information may be incomplete or unable to be used in any situation. The information in this report should not be interpreted by artificial intelligence software. No reliance or actions must therefore be made on that information without seeking prior expert professional, scientific and technical advice. To the extent permitted by law, Genics (including its affiliates, directors, employees and consultants) excludes all liability to any person for any consequences, including all losses, damages, costs, expenses and any other compensation, arising directly or indirectly from using this publication (in part or in whole) and any information or material contained in it.



Genics Oral Fluid Collection Guide



GENICS

genics.com

Hello

Welcome to **Genics Oral Fluid Collection Guide**. Here you will find a list of all materials needed to complete the collection including clear instructions on the process.

You will need

- Cotton rope (approx. 100 cm long)
- Plastic screw top tube 50ml
- Clean plastic bags
- Disposable gloves
- Permanent marker
- Submission form
- Box cutter, knife or scissors



Image 1) With a rope (approx. 70cm long) tie a half hitch knot or similar then unravel the strands.

Steps

- 1 Sample every 3rd pen per barn as shown below. Take samples at times when animals are most active

1	2	3	4	5	6	7	8	9
18	17	16	15	14	13	12	11	10

- 2 Wearing the gloves provided, tie a half hitch knot (see Image 1 left) securing the rope in a clean area of the pen away from feed and water. To avoid contamination, it is important the rope does not come in contact with the floor



Image 2) Strip the rope so the fluid accumulates in the corner of the bag

- 3 For pens with 25 or more pigs, two or more ropes may be used for collection (samples may be pooled by pen)
- 4 Ropes should remain in pens for 20-30 minutes. While waiting, label tubes consecutively (e.g: Tube #3,6,9 etc)
- 5 Wearing gloves, collect the ropes and insert the wet end into a clean plastic bag
- 6 Strip the rope so fluid accumulates in the corner of the bag (See Image 2 middle). To avoid cross contamination, change gloves between samples
- 7 Cut a corner of the plastic bag and drain fluid in the 50ml tube provided
- 8 A minimum sample of 2.5ml and maximum of 40ml is required
- 9 Discard all ropes after collection. Never reuse plastic bags or collection tubes
- 10 Do not pool oral fluid samples from ropes in different pens.

Shipping instructions overleaf

